



A Fast, Compact TPC with GEM Readout for Tracking and Electron Identification at RHIC

Craig Woody

Brookhaven National Lab

Linear Collider Workshop
Berkeley, CA

October 18, 2003

Participants

BNL

- Physics

C.Woody, B.Azmoun, M.Sivertz, J.Mitchell, M.Purschke

- Instrumentation

P.O'Connor, B.Yu, V.Radeka, G.Smith, J-F.Pratte

Yale

N.Smirnov, R.Majka

Nevis

C-Y.Chi, W.Sippach

Tokyo

K.Ozawa, T.Isobe, M.Inuzuka, H.Hamagaki, T. Sakaguchi,
F.Kajihara, T.Gunji, N.Kurihara, S.Sawada (KEK), S.Yokaichi (RIKEN)

Weizmann




I.Tserruya, I.Ravinovich, L. Shekhtman, A.Kozlov, Z.Fraenkel

Stony Brook

A.Milov



PHENIX and STAR at RHIC



PHENIX and STAR Upgrade Plans

Extended Physics goals

- ❑ Particle Identification at high Pt
- ❑ Charm/Beauty measurements
- ❑ Low Mass Dilepton Pairs
- ❑ High Pt and Jets

PHENIX Detector upgrades

- Aerogel Cerenkov detector
- Silicon Vertex Tracker
- **HBD and inner TPC**
- Enhanced muon trigger


STAR Detector Upgrades

- DAQ upgrade ($\rightarrow 100 \text{ Hz} \rightarrow \text{KHz}$)
- Upgrade readout electronics
- RPC TOF detector
- Silicon Active Pixel Vertex detector
- **New fast, compact TPC**

Physics measurements addressed by the Fast, Compact TPC in PHENIX

- Low mass lepton pairs and vector mesons
- Charm and B physics with resolved secondary vertices
 - low mass tracking just outside vertex detector
 - allows measurement in both heavy ion and pp running
- Improved inner tracking for PHENIX
 - increased h and f coverage
(needed for jet and g-jet physics)
 - tracking through the magnetic field
(improves momentum resolution, the ability to measure real low p_T tracks, and to reject high p_T background tracks)


TPC/HBD Detector




GEMs are used for both TPC and HBD

- Fast, compact TPC
 $R < 70 \text{ cm}$, $L < 80 \text{ cm}$, $T_{\text{drift}} \leq 4 \text{ msec}$
- Serves as an inner tracking detector in both HI and pp, providing tracking through the central magnetic field
 $Df = 2p$, $|h| \leq 1.0$
 $Dp/p \sim .02p$
- Provides electron id by dE/dx
 $\Rightarrow e/p$ separation below 200 MeV
- HBD is a proximity focused Cherenkov detector with a $\sim 50 \text{ cm}$ radiator length
- Provides electron id with minimal signals for charged particles
 \Rightarrow "Hadron Blind Detector"


PHENIX Inner Magnetic Field



\pm Configuration




Inner field itself
is non-uniform



Inner Coil
creates a
“field free”
($B_{dl}=0$)
region inside
the Central
Magnet


Field Integrals



Hadron Blind Detector

Multistage GEM detector with CsI Photocathode

- CsI photocathode deposited on outer GEM foil (must first deposit Ni+Au on GEM foil)
- Multistage GEM used to detect few tens of photoelectrons
- Minimal ionization for charged particles in region above first GEM foil
- Fairly coarse segmentation of readout plane ($\sim 10 \text{ cm}^2$ pads)



Choice of Gases

HBD

- Need deep UV transparency for highest photoelectron yield
- Windowless design \Rightarrow same gas for radiator as working gas for GEMs

Gas	Ecutoff (eV)	N_0 (cm^{-1})	dE/dx (keV/cm)
CF ₄	11.5	936	7
CH ₄	8.5	185	1.5
Ne	15	2664	1.6
Ar	9	255	2.4

TPC

- Fast drift velocity (~ 10 cm/msec)
- Low diffusion
- High dE/dx
- Same drift gas as working gas for GEMs


CF₄ (or mixtures
containing CF₄) look
most promising

Would like to find a single gas (or gas mixture) which would work as a drift gas, radiator gas and working gas for GEMs


HBD R&D at Weizmann




See talk by I.Ravinovich
(N25-3, Gas Detectors II, Wed, 10/22, 10:30-12:00)



Triple GEM operating in CF₄
with CsI photocathode




PHENIX TPC Tracking Simulation



Currently
developing TPC
tracking and
pattern
recognition
software using
existing PHENIX
tracking
detectors and
future VTX
detector

Event display of the PHENIX TPC for a
central Au-Au event
(PISA simulation - J.Mitchell)

TPC Readout




Readout parameters

Number of pads	80K
Pad size	2x10 mm ²
Drift time	3.5 msec
Sampling rate	40 MHz (25 ns)
Sampling resolution	2.5 mm, 8 bits
Number of samples	140
Unsuppressed data volume	11 MB
Suppressed data volume (~1/10)	1 MB
Actual data volume	100 KB
Buffer latency	4 msec
Readout time	40 msec
Data transmission rate	200 Gbit/sec
Power per channel	100 mW
Total power	8 KW

Requirements

- Very high density readout
- Need to minimize power
- Low noise, zero suppression
- Minimize data volume (triggering)

TPC/HBD readout electronics




- High density, low power preamp and shaper (custom ASIC)
- 40-50 MHz, 8-10 bit FADC (TI, AD - Nevis)
- Digital baseline subtraction, zero suppression
- Provide buffering for Level 1 trigger (~ 5 events)

Similar 32 ch preamp/shaper being developed for medical imaging
0.18 mm CMOS, ~ 4 mW per channel
Final ASIC 4.3 x 1.6 mm²

See talk by J-F. Pratte (N1-7, Analog & Digital Circuits I, Monday, 10/20, 13:30-15:15)

Interpolating Pad Readout



Two Intermediate Strips



Single Intermediate Zigzag

B. Yu et.al., IEEE Trans. Nucl. Sci.
Vol. 50, No. 4. (2003) 836-841




Fine “Zigzag” pattern



Overall position error: 93μm rms

Including ~ 100μm fwhm x-ray p.e. range,
100μm beam width, alignment errors

Test Drift Cell




- Used to study**
- Drift velocities
 - Drift lengths
 - Diffusion parameters
 - Energy loss (dE/dx)
 - Study impurities
 - Readout structures
 - Field cage design

Lab 2-86 in Physics

Drift Cell Measurements


Ar/CO₂ (80/20)





22% FWHM/mean

Energy Resolution
(Top ⁵⁵Fe source)

Ar/CF₄ (95/5)



28% FWHM/mean




Drift Velocities (laser)


Charge Distribution from Drift Cell

(top source, 2 mm pad readout)

Ar/CO₂ (80/20)




Ar/CF₄ (95/5)




Measured

S=1.1 mm





Simulation
(N.Smirnov)

S=1.4 mm



Testing Small 3M Foils

(courtesy of J. Collar, Univ. of Chicago)




80 **mm** (55 **mm**) holes spaced in a hexagonal pattern with 140 **mm** pitch (approx. same as CERN foils)

Visually looks to be excellent quality


3 stage GEM with three 1" dia. 3M foils
(2 x 1.5 mm transfer gaps + 2 mm induction gap)

Comparison of Small Triple GEMS

Ar/CO₂ (80/20)





CF₄




3M GEMs give higher gain and have somewhat different slopes

Gain Stability - Small 3M GEMs

Ar/CO₂ (80/20)



CF₄




We observe a large gain variation with time with the small 3M GEMs

Time scales involved imply a possible charging effect in the foils


(Gas: 35 ppm H₂O, <1ppm O₂)

Gain Stability - Small CERN GEMs in Ar/CO₂




No large gain variations observed with small CERN GEMs

Large (10x10 cm²) 3M GEM Foils




- 1st batch - “undercut” (faster etching kapton \Rightarrow “over etched”)
- 2nd batch - “conical” (standard kapton but smaller holes)
(70mm \rightarrow 40 mm vs 85 mm \rightarrow 55 mm)


Comparison of Gain Curves for 10x10 cm² Foils



Gain Stability - Large Foils



TPC Prototype - To be built this year



- Full size module
- Tests field cage design
- Tests materials to be used in actual construction
- Allows testing with cosmic rays
- Provides structure on which to design and study readout plane
- Test bench for new readout electronics

Summary

- We have an ambitious goal to build a single detector that will serve many purposes:
 - Fast drift, high resolution TPC
 - Hadron Blind Cherenkov Detector
 - Operate in high multiplicity heavy ion and high luminosity pp collisions
- Requires meeting new challenges in detector technologies
 - GEM detectors w/CsI photocathodes
 - Operation with CF_4
 - New, high density, low power readout electronics
 - High data volumes and data rates
- We need a supplier of good quality GEM foils for all future development work, and are interested in working with all potential vendors